

# DIGITAL BATHROOM SCALE

## INSTRUCTION MANUAL

### Key Features

- High precision strain gauge sensors system for accuracy
- Large LCD display with backlight
- Tempered black glass platform
- Three different units of measure: lb/kg/st
- Capacity: 400lbs / 200kg
- Auto-on, Auto-zero, Auto-off
- Low battery and Overload indication
- Powered by 2 x AAA battery (included)

### Attention

- Remove shoes before stepping on the scale to avoid scraping or scratching the glass platform.
- Warning: the glass platform may be slippery when wet, so be sure both the glass and your feet are dry before using the scale.
- Use soft tissue with glass cleaner to clean the surface of the scale. Do not use soap or other chemicals.
- Keep the scale away from water, corrosive liquid, heat and extreme coldness.
- The scale is a high precision measuring device. Please handle with care to avoid dropping or breaking.
- The scale is not suitable for professional use.
- Please check the battery if the scale malfunctions. Change the battery when "LO" appears on the display which means the battery is low.

- "Err" showed on display indicates overload.



### Setting up your scale

- Insert 2 AAA battery into the battery compartment on the back of the scale and replace the cover.
- Press the button on the back of the scale to select measure units (lb/kg/st).

### How to use the scale

- Gently press on the scale with your foot to turn it on. Display will show "0.0".
- Step on the scale and keep your body still until the scale stabilizes your weight.
- The scale will indicate your weight in the unit of measure selected (lb/kg/st).
- Step off the scale. The scale will automatically turn off.

